

DOCTOR _____ ACCT# _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

PHONE _____ EMAIL _____

PATIENT NAME _____

DATE SHIPPED _____ DATE NEEDED* _____

**Date needed should be at least 1 day before appointment date.*

APPROVAL TO CHARGE EXPEDITED MANUFACTURING & SHIPPING

DIGITAL SCAN TAKEN WITH:

iTero® Carestream CEREC


TRIOS® Medit Other _____

OFFICE USE: 1 2 3 4 + PD: SA DR

MODELS: U L BOTH BANDS CROWNS BROKEN

IMPRESSIONS: U L BOTH

DISINFECT: _____ QA IN: _____ FINAL INSP: _____



CLEAR IMAGE® ALIGNERS

Upper Lower

Add Digital Pontics

Refinement

Upper Lower

Location _____

Shade (Paint) _____

Send Tx Review*

Yes* Email _____

Opposing Model Required for Warranty

*By selecting this you will be responsible for a planning fee if patient does not proceed with treatment.

CHOOSE TEETH TO RESET

R	7	6	5	4	3	2	1	1	2	3	4	5	6	7	L
	7	6	5	4	3	2	1	1	2	3	4	5	6	7	

Reset all teeth as needed Reset anterior only Reset 5-5 only

SPACE & OVERJET INSTRUCTIONS

- Okay to leave space on upper due to overjet restrictions
Leave space _____
- Opposing arch unavailable. Overjet measurement is: _____

ATTACHMENT INSTRUCTIONS

- No attachments
- Attachments only on: _____
- Size/Shape preferred: _____
- Remove attachments on last stage

IPR INSTRUCTIONS

- Reduce Teeth as Needed in Laboratory
- No IPR Required

OTHER INSTRUCTIONS

- I want to correct _____
- Midline _____
- Overbite _____
- Overjet _____
- Spacing _____

PONTIC OPTIONS

- Add Pontics
- Location _____ Shade _____

INVISIBLE RETAINERS

	U	L
Single Invisible Retainer	<input type="radio"/>	<input type="radio"/>
Guardian™ 2	<input type="radio"/>	<input type="radio"/>
Guardian™ 3	<input type="radio"/>	<input type="radio"/>
Guardian™ 4	<input type="radio"/>	<input type="radio"/>
IR Express <small>*5 Day Turnaround. See exclusions on page 2</small>	<input type="radio"/>	<input type="radio"/>

FIXED LINGUAL RETAINERS (FLR)

	U	L
Placement of Retainer		
Central - Central	<input type="radio"/>	<input type="radio"/>
Lateral- Lateral	<input type="radio"/>	<input type="radio"/>
Cuspid - Cuspid	<input type="radio"/>	<input type="radio"/>
Placement of Pads		
Composite Pads on Each Tooth	<input type="radio"/>	<input type="radio"/>
Composite Pads on Distal Most Teeth	<input type="radio"/>	<input type="radio"/>
Mesh Pads on Each Tooth	<input type="radio"/>	<input type="radio"/>
Mesh Pads on Distal Most Teeth	<input type="radio"/>	<input type="radio"/>
Type of Wire		
Round .028	<input type="radio"/>	<input type="radio"/>
.016 x .022 Braided	<input type="radio"/>	<input type="radio"/>
.016 x .022 Solid Stainless Steel	<input type="radio"/>	<input type="radio"/>

SPECIAL INSTRUCTIONS (For aligners, please describe treatment goals and patient chief complaint in detail)

DOCTOR SIGNATURE _____

License # _____ Expiration _____

(800) 522-4636
www.specialtyappliances.com

CLEAR IMAGE® ALIGNERS

Clear Image® Aligners are processed from your scans, models, or impressions. Aligner specialists will determine the optimal number of aligners required for treatment if you do not specify an amount. Aligners are made using Zendura™ material. This product can be ordered in single or dual arch.

Tx Review

All cases require doctor approval of the setup prior to aligner fabrication. If the setup is not approved or the patient does not proceed with treatment, only the Treatment Planning Fee and any applicable model and shipping fees will be incurred.

Treatment Planning Fee

The Treatment Planning Fee is waived for single arch cases with 10 or more aligners and dual arch cases with 20 or more aligners.

Choose Teeth to Reset

The diagram provided allows the clinician to select individual or groups of teeth for corrections by marking the diagram accordingly.

Opposing Models

Specialty Appliances requires an opposing scan or model be submitted when making Clear Image® Aligners. This ensures no occlusion interference is created by tooth movements and guides us on necessary tooth movements for best outcomes. If an opposing scan or model is not submitted, the warranty for success of Clear Image Aligners is not offered and any refinement aligners will be subject to full fees.

IPR Instructions

Most Clear Image® Aligners require IPR to help create space for the teeth to move. In this section, please specify whether Specialty Appliances is to determine IPR amounts and locations or that no IPR is required.

Scanning Stone Models

Clear Image® Aligner cases that are sent in as impressions or stone models will be subject to automatic impression pouring and scanning fees. If the patient chooses not to proceed with treatment, these fees will still apply.

GUARDIAN® INVISIBLE RETAINERS

GUARDIAN® Retainers are processed from your scans, models, or impressions. From your records, we will produce a printed model for laboratory fabrication at no charge. GUARDIAN® Retainers are available in quantities 2,3, and 4 per arch and are fabricated all at one time. All GUARDIAN® Retainers are processed on printed plastic models for maximum accuracy. These models are returned to give your practice the option to fabricate replacement retainers as needed. Zendura™ is used as a standard for all GUARDIAN® Retainers.

GUARDIAN® ENHANCEMENTS

Remove Brackets

We routinely remove all brackets on GUARDIAN® cases unless requested otherwise. Using the latest in advanced software, all fixed appliances are removed, free of charge, at the scanning stage to be reflected in the printed model. We recommend that all wires and bands are removed prior to submitting scans or impressions. If wires or bands are present, it will result in additional charges to remove them.

Add Fixed Lingual Retainer to a GUARDIAN® Series

Specialty Appliances will build a Fixed Lingual Retainer (FLR) design of your choice and process the GUARDIAN® Retainers to fit over the FLR. At the appointment when the braces are removed, the FLR is delivered with the GUARDIAN® Retainers fitting precisely over the FLR. Several designs of Fixed Lingual Retainers are available. Please use the diagram to indicate your wire and pad preferences.

Add Pontics to GUARDIAN® Retainers

Pontics can be added to any GUARDIAN® Retainer. We can also add pontics to any FLR design as a semi-permanent solution for patients too young for implants or restorations.

We encourage practices to contact our customer service representatives for additional information regarding these options.

IR EXPRESS

Bracket removal and shipping costs are included with IR Express. IR Express is fabricated with Zendura A .030 and has a 5 day turnaround time. IR Express does not include a printed model.

Exclusions:

Cases that require band, wire or expander removal, pontics, bonded retainers, anterior or posterior bite planes or Theroux do not qualify for IR Express. Material other than Zendura A .030 does not qualify for IR Express. Cases that are submitted analog (stone models/impressions) do not qualify for IR Express.