

Ortho SurgiArch™

SurgiArch™ is a patent-pending 3D-printed wire bonded to the teeth prior to oral surgery to set arch position.

SurgiArch Features

- ◆ Cross section is created at each tooth for retention for bonding
- ◆ Seating rests to help seat surgical wire in proper position
- ◆ Hooks allow surgeon to tie arches together for healing
- ◆ 4 designs for pads at cross section

SurgiArch Uses

- ◆ Trauma from fall or automobile accidents
- ◆ Sleep Apnea surgeries
- ◆ Orthognathic Surgery Before Aligner Treatment



X-Bar



Cross Bar



Mesh Pad



Perforated Mesh Pad